

# MODULAR MEDICATION CASSETTE

**DRUG TRANSFER  
STORAGE SYSTEM**

**EGASUS**  
MEDICAL CONCEPTS INC.





## UNIQUE FEATURES

- Multiple configurations
- Easy to transport
- Up to 3 levels with a maximum of 18 bins
- Constructed of lightweight plastic (ABS)

## SAFETY

- Various security lock solutions:
  - Key lock
  - Electronic security code

## THE STANDARD IN MODULAR STORAGE SYSTEMS

- A Compatible with ISO-3394
- 60 x 40cm

## FLEXIBILITY

- 2 different orientations:
  - To distribute medication: 60cm front
  - To transport in trolleys: 40cm front
- Bins of two interchangeable sizes

## DIN 55 510 / ISO3394 / HTM71

The Pegasus medical modular medication cassette has a footprint of 60 X 40cm, therefore, it can easily slide into any cart and cabinet that are designed to accommodate ISO 3394 modular systems.

The Pegasus medical modular medication cassette can slide into a side panel of an ISO 3394 cart or cabinet in 2 different orientations: 60cm front for medication distributions purpose & 40cm front for transfer carts applications.



## EXAMPLE OF A LOGISTIC CHAIN WITH THE USE OF A MEDICATION CASSETTE



The modular medication cassette is perfect for use in a drug distribution process

← Used cassette will be transferred back to the pharmacy

→ Newly loaded cassettes from the central pharmacy are distributed to the departments where the cassettes are being exchanged full to empty.



### ELECTRONIC CODE LOCK

A user code can easily be programmed after the entry of a master code



### KEY LOCK

Per cassette, we deliver 1 set of 2 keys





### INFORMATION

#### Cleaning

The cassettes, medication drawers and locks can be cleaned by hand with PH neutral cleansers

#### Delivery



Keys are always delivered in one set of two keys

<b>Cassettes</b>				
<b>Item Number</b>	<b>MED-CAST-1L</b>	<b>MED-CAST-2L</b>	<b>MED-CAST-3L</b>	<b>MED-CAST-4L</b>
<b>Dimensions</b>	14.8x59.5x40cm	24.3x59.5x40cm	33.9x59.5x40cm	43x60x40cm
<b>Material</b>	ABS	ABS	ABS	ABS
<b>Drawers Configuration</b>	4 large drawers or 6 small drawers	8 large drawers or 12 small drawers	12 large drawers or 18 small drawers	16 large drawers or 24 small drawers

## SELECT THE LOCK TYPE AND ADD TO THE END OF THE ITEM NUMBER

<b>Locking System</b>			
	<b>Without</b>	<b>Key lock</b>	<b>Electronic code lock</b>
<b>Item Number</b>	<b>...-W</b>	<b>...-WKL-W</b>	<b>...-WCL-W</b>

<b>Medication drawer</b>				
<b>Item Number</b>	<b>MED-DWR-80 (white)</b>	<b>MED-DWR-120 (white)</b>	<b>MED-DWR-120C (clear)</b>	<b>MED-DWR-80C (clear)</b>
<b>Dimensions</b>	8x34x8cm	12x34x8cm	12x34x8cm	8x34x8cm
<b>Material</b>	ABS	ABS	ABS	ABS

<b>Medication drawer dividers</b>		
<b>Item Number</b>	<b>MED-DWRDIV-80</b>	<b>MED-DWRDIV-120</b>
<b>Dimensions</b>	7x5cm	11x5cm
<b>Material</b>	ABS	ABS